

## LEICA CM3050 S CRYOSTAT



<b>Specimen</b>	<ul style="list-style-type: none"> <li>Maximum specimen size: 40x55mm</li> <li>Overall specimen feed : 25mm</li> </ul>
<b>Temperatures</b>	<ul style="list-style-type: none"> <li>Specimen temperature:-10°C/-50°C</li> <li>Cryochamber temperature: 0/-40°C</li> </ul>
<b>Section thickness</b>	<ul style="list-style-type: none"> <li>Sectioning: 0.5-300µm <ul style="list-style-type: none"> <li>0.5-2 µm (in 0.5 µm steps)</li> <li>2-10 µm (in 1 µm steps)</li> <li>10-20 µm (in 2 µm steps)</li> <li>20-50 µm (in 5 µm steps)</li> <li>50-100 µm (in 10 µm steps)</li> <li>100-300 µm (in 50 µm steps)</li> </ul> </li> <li>Trimming: 5, 10, 30, 50, 100, 300 µm</li> </ul>
<b>Sectioning options</b>	<ul style="list-style-type: none"> <li>Manual</li> <li>Motorized: Intermittent stroke; Continuous stroke; Single stroke</li> </ul>
<b>Knife holder and blades</b>	<ul style="list-style-type: none"> <li>Leica CE with anti-roll guide</li> <li>Disposable low profile blades</li> </ul>
<b>Accessories</b>	<ul style="list-style-type: none"> <li>Footswitch     Stationary heat extractor</li> <li>90° Prism for low temperature (-50°C) quick freezing</li> <li>Thermal plate Leica (HI1220)</li> </ul>