

LEICA TP1020 AUTOMATIC TISSUE PROCESSOR *



Tissue basket (1)	<ul style="list-style-type: none"> Capacity: <ul style="list-style-type: none"> 75 cassettes (22x23x4mm) 38 Mega cassettes (27x23x10mm)
Reagent stations (9)	<ul style="list-style-type: none"> Capacity: 1.8L
Paraffin stations (2)	<ul style="list-style-type: none"> Vogel Histo-comp paraffin wax 56°C Capacity: 1.8L Temperature setting range: 45-65°C
Operation:	<ul style="list-style-type: none"> Manual Automatic: <ul style="list-style-type: none"> 9 programs, freely selectable Delayed start time of processing (up to 9 days)
Accessories	<ul style="list-style-type: none"> Fume cupboard and material for preparation of reagents and samples Disposal of waste reagents according to UAB normative Leica EG 1150H Paraffin Embedding Station Leica EG1150C Cold plate

Temporally donated by Professor Juan Hidalgo's group, only to be used by the histology lab technician