

## SME CRYOTOME Thermo Electron Corporation



<b>Specimen</b>	<ul style="list-style-type: none"> <li>• Maximum specimen size: 40x55mm</li> <li>• Overall specimen feed : 25mm</li> </ul>
<b>Temperatures</b>	<ul style="list-style-type: none"> <li>• Specimen temperature: 0°C/-40°C</li> <li>• Cryochamber temperature: 0/-50°C</li> <li>• Cryobar: -20/-40°C</li> </ul>
<b>Section thickness</b>	<ul style="list-style-type: none"> <li>• Sectioning: 0.5-60 µm <ul style="list-style-type: none"> <li>○ 0.5-20 µm (in 1 µm steps)</li> <li>○ 20-60 µm (in 5 µm steps)</li> </ul> </li> </ul>
<b>Sectioning options</b>	<ul style="list-style-type: none"> <li>• Manual</li> <li>• Motorized: <ul style="list-style-type: none"> <li>○ Continuous stroke</li> <li>○ Single stroke</li> </ul> </li> </ul>
<b>Knife holder and blades</b>	<ul style="list-style-type: none"> <li>• 7408 with anti-roll guide</li> <li>• Disposable low profile blades</li> </ul>
<b>Accessories</b>	<ul style="list-style-type: none"> <li>• Footswitch</li> <li>• Stationary heat extractor</li> <li>• Fumigation</li> <li>• Thermal plate Leica (HI1220)</li> </ul>