


DAKO AUTOSTAINER PLUS



Samples	<ul style="list-style-type: none"> Up to 48 microscope slides (25 x 75 x 1mm) <ul style="list-style-type: none"> Paraffin embedded tissue sections Frozen tissue sections Cytospins Cell smears
Reagents	<ul style="list-style-type: none"> Up to 64 different reagents per run (15ml/vial) Ten liter buffer container Dispense volumes: 100, 150, 200, 400 and 600 µl
Software options	<ul style="list-style-type: none"> Different protocols per staining run Up to 35 steps in each protocol (including rinse steps) Incubation time: 2-90 minutes Dispense volume (for each slide and step) Reagent dispense location (for each slide)
Accessories	<ul style="list-style-type: none"> Seymour glass label printer HP inkjet printer Fume cupboard and material for preparation of reagents and samples Disposal of waste reagents according to UAB normative