

CONFOCAL LASER SCANNING MICROSCOPE ZEISS LSM 700:



4 Laser lines	<ul style="list-style-type: none">• Laser line 405 nm (5mW fiber output)• Laser line 488 nm (10mW fiber output)• Laser line 555 nm (10mW fiber output)• Laser line 639 nm (5mW fiber output)
Objectives	<ul style="list-style-type: none">• W N-Achroplan 10x/0,3 M27 objective (dt=2,6 mm)• Plan-Apochromat 20x/0.8 Ph2 M27 objective (dt=0.55mm)• W Plan-Apochromat 40x/1,0 Ph3 M27 objective (dt=2,5mm), VIS-IR• Plan-Apochromat 40x/1,3 Oil DIC M27 objective (dt=0,2 mm), (UV)VIS-IR• Plan-Apochromat 63x/1,40 Oil DIC M27 objective (dt=0,19 mm)
Maximum lateral resolution	<ul style="list-style-type: none">• 200 nm
Detection system	<ul style="list-style-type: none">• Two confocal fluorescence detectors (High sensitivity photomultipliers)• Software ZEN
Applications	<ul style="list-style-type: none">• 3D Image• Multiple fluorescence and colocalization analysis• <i>In vivo</i> image acquisition (time-lapse)

OPTICAL EPIFLUORESCENCE MICROSCOPE NIKON ECLIPSE 90i:



Illumination	<ul style="list-style-type: none">• Visible DC lamp 12V, 100W• Mercury lamp 100W
Objectives	<ul style="list-style-type: none">• Objective Plan Apo 2x/0.1• Objective Plan Fluor 10x/0.30 DIC L/N1• Objective Plan Apo 20x/0.75 DIC M/N2• Objective Plan Apo 40x/0.95 DIC M/N2• Objective Plan Apo 60x/1.4• Objective Plan Apo 100x/1.4
Magnification	<ul style="list-style-type: none">• 16-2000x
Fluorescence filters	<ul style="list-style-type: none">• DAPI: Ex 340-380 / DM 400 / BA 435-485• FITC: Ex 465-495 / DM 505 / BA 515-555• G-2A: Ex 510-560 / DM 575 / BA 590
Detection system	<ul style="list-style-type: none">• Nikon Camera DXM1200F (Software ACT-1)
Observation methods	<ul style="list-style-type: none">• Brightfield• Epifluorescence• DIC or Nomarski• Polarized light

OPTICAL EPIFLUORESCENCE MICROSCOPE NIKON ECLIPSE 80i:



Illumination	<ul style="list-style-type: none">• Visible DC lamp 12V, 100W• Mercury lamp 100W
Objectives	<ul style="list-style-type: none">• Objective Plan UW 2x/0.06• Objective Plan Fluor 4x/0.13• Objective Plan Fluor 10x/0.30 DIC L/N1• Objective Plan Apo 20x/0.50 DIC M• Objective Plan Apo 40x/0.95 DIC M/N2• Objective Plan Apo 100x/1.4 Oil DIC H
Magnification	<ul style="list-style-type: none">• 20-1000X
Fluorescence filters	<ul style="list-style-type: none">• UV-2A: Ex 330-380 / DM 400 / BA 420• B-2A: Ex 450-490 / DM 505 / BA 520• G-2A: Ex 510-560 / DM 575 / BA 590
Detection system	<ul style="list-style-type: none">• Nikon Camera DXM1200F (Software ACT-1)
Observation methods	<ul style="list-style-type: none">• Brightfield• Epifluorescence

OPTICAL EPIFLUORESCENCE MICROSCOPE NIKON ECLIPSE TE2000-E:



Illumination	<ul style="list-style-type: none">• Visible DC lamp 12V, 100W• Mercury lamp 100W
Objectives	<ul style="list-style-type: none">• Objective Plan Fluor 10x/0.30 Ph1 DLL• Objective Plan Apo 20x/0.45 DIC Ph1 DM• Objective Plan Apo 40x/0.60 DIC M• Objective S Fluor 40x/1.3 Oil• Objective Plan Apo 100x/1.3 Oil Ph3 DLL
Magnification	<ul style="list-style-type: none">• 100-1500X
Fluorescence filters	<ul style="list-style-type: none">• DAPI: Ex 340-380 / DM 400 / BA 435-485• FITC: Ex 465-495 / DM 505 / BA 515-555• G-2A: Ex 528-553 / DM 565 / BA 578-633
Detection system	<ul style="list-style-type: none">• Nikon Color Camera DS-2Mv (Software NIS-Elements)• Hamamatsu Monochromatic Camera ORCA-ER (Software Metamorph)
Observation methods and applications	<ul style="list-style-type: none">• Brightfield• Epifluorescence• Phase contrast (10x and 20x)• Fixed or <i>in vivo</i> samples, with temperature controller for RT 50°C• <i>Stacks</i> in z and <i>time-lapse</i> (motorized stage)

OPTICAL EPIFLUORESCENCE MICROSCOPE NIKON ECLIPSE TE2000-U:



Illumination	<ul style="list-style-type: none"> • Visible DC lamp 12V / 100W • Cairns monochromator with Xenon lamp 150W, wavelengths between 300-700 nm
Objectives	<ul style="list-style-type: none"> • Objective Plan Fluor 10x/0.30 Ph1 DLL • Objective Plan Apo 20x/0.45 DIC Ph1 DM • Objective Plan Apo 40x/0.60 DIC M • Objective S Fluor 40x/1.3 Oil • Objective S Fluor 100x/0.5-1.3 Oil
Magnification	<ul style="list-style-type: none"> • 100-1500X
Fluorescence filters	<ul style="list-style-type: none"> • D535/25 (400-600 nm) • D605/55M (500-700 nm) • D510/40M (320-620 nm) • HQ530/50 (300-750 nm)
Detection system	<ul style="list-style-type: none"> • Hamamatsu Camera ORCA-ER (Software MetaFluor)
Observation methods and applications	<ul style="list-style-type: none"> • Brightfield • Epifluorescence • Phase Contrast (10x and 20x) • <i>In vivo</i>, with temperature controller for RT 50°C • <i>Time-lapse</i>: intracellular calcium, pH, membrane potential.....