




CONFOCAL LASER SCANNING MICROSCOPE ZEISS LSM 700:



| 4 Laser lines | Laser line 405 nm (5mW fiber output) |
|----------------------------|--|
| | Laser line 488 nm (10mW fiber output) |
| | Laser line 555 nm (10mW fiber output) |
| | Laser line 639 nm (5mW fiber output) |
| Objectives | W N-Achroplan 10x/0,3 M27 objective (dt=2,6 mm) |
| | Plan-Apochromat 20x/0.8 Ph2 M27 objective |
| | (dt=0.55mm) |
| | W Plan-Apochromat 40x/1,0 Ph3 M27 objective |
| | (dt=2,5mm), VIS-IR |
| | Plan-Apochromat 40x/1,3 Oil DIC M27 objective |
| | (dt=0,2 mm), (UV)VIS-IR |
| | Plan-Apochromat 63x/1,40 Oil DIC M27 objective |
| | (dt=0,19 mm) |
| Maximum lateral resolution | • 200 nm |
| | Two confocal fluorescence detectors |
| Detection and acquisition | (High sensitivity photomultipliers) |
| system | Software ZEN |
| | |
| Applications | 3D Image |
| | Multiple fluorescence and colocalization analysis |
| | In vivo image acquisition (time-lapse) |
| 1 | |